


Omnipod DASH[®] System Guide for Automatic Data Uploads

Track trends in your diabetes care and share your Omnipod DASH[®] System data with your care team.




Connect PDM to Wi-Fi

Data automatically
uploaded once every 24 hours

Data available in Glooko
for you and your care team
to view and interpret

A step-by-step guide to:

- Create a PodderCentral[™] Account
- Create a Glooko Account
- Connect to a secure Wi-Fi network
- Enable automatic upload and share data



To prepare for initial set up, make sure you have access to:

- a computer or smartphone
- your PDM serial number
- a secure WiFi network
- personal email inbox


I. Create a PodderCentral™ Account

- 1 Go to **PodderCentral.com**.
 - a. On your computer or smartphone web browser, enter the email address you provided when you placed your first Omnipod® order.

NOTE: Internet Explorer is not supported. Use Chrome, Safari or FireFox. Avoid having multiple browsers open while creating your account.
 - b. You may be presented with an option to select **“Product Source”**.
 - If you received your supplies from anyone outside of Kaiser, please select **“Other”**.
 - c. You will be presented with a ‘Thank You’ message directing you to check your inbox for your invitation to sign up for PodderCentral™.
 - d. Select **“Close”** and you will be returned to the PodderCentral™ login page.
- 2 You will now need to leave the PodderCentral™ page and go to the inbox of the email you provided.
 - a. Locate the PodderCentral™ email invitation with the subject line **“You’re Invited to PodderCentral™”**.
 - b. Select the **“unique link”** within the email.
- 3 You will be brought to the PodderCentral™ **“Account Setup 1 of 2”** and must fill in the requested user information.

TIP: When inputting information, make sure to type in individual fields. Do not auto-fill.

 - a. Once all information has been completed, select **“I’m not a robot”**.
- 4 You will be brought to **“Account Set up 2 of 2”**. You must answer three security questions from the drop-down menu and certify that you:
 - are 18 years of age or older;
 - have read and understood Insulet Corporation’s HIPAA Privacy Notice, Privacy Policy, and Terms of Use;
 - authorize Insulet Corporation, its distributors, affiliates and wholly-owned subsidiaries to contact you by phone or e-mail regarding the Omnipod® Insulin Management System and other diabetes-related supplies and services.
- 5 Select **“Register”**.



PodderCentral™ | glooko

Log into PodderCentral™ to Automatically Upload DASH™ Data to Glooko

* Required Fields

Username *

Forgot Username ☐ Remember Username

Password *


Forgot Password

LOG IN

Or Create a PodderCentral™ Account

You will need a PodderCentral™ account to share your data with Glooko. Please enter your Email address to begin sign up.

Step 1



PodderCentral™ | glooko

Hello Podder™,

We'd like to invite you to join PodderCentral™, an online portal just for our Podders™ that puts control in your hands. Once you activate your PodderCentral™ account, you're just a few clicks away from being able to:


- + Automatically upload your data to Glooko

Your account information is now at your fingertips, where and when you want it. Access your account using the My Omnipod® Mobile App or your computer, it's up to you!

Use this [unique link](#) now to get started and activate your PodderCentral™ account.

Sincerely,

Step 2



PodderCentral™ | glooko

* Required Fields

Account Setup

1 of 2

Please enter information for the person wearing the Pod

Omnipod® User First Name *

Omnipod® User Last Name *

Omnipod® User D.O.B.*


MM DD YYYY

Username *

Password *

Product Source

Step 3



PodderCentral™ | glooko

* Required Fields

Account Setup

2 of 2

Security Question 1 *

Answer *

Security Question 2 *

Answer *

Security Question 3 *

Answer *

I, clara.hill certify that:

☐ I am 18 years of age or older

☐ I have read and understand Insulet Corporation's HIPAA Privacy Notice, Privacy Policy and Terms of Use

☐ I authorize Insulet Corporation, its distributors, affiliates and wholly-owned subsidiaries to contact me by telephone or e-mail regarding the Omnipod Insulin Management System and other diabetes-related supplies and services.

Step 4

II. Create a Glooko Account*

- 1 Open a web browser and go to join.glooko.com.
 - a. Enter your email address; password; first name and last name of patient; country of residence; birthdate of patient; and phone number (*optional*).

TIP: Make sure to add **"US"** for country when setting up your Glooko account.
 - b. Select if you'd like to receive text messages from Glooko.
 - c. Enter your Healthcare Provider's ProConnect Code.

TIP: Refer to Section VI in this guide for more information on how to connect the ProConnect code to your account.

NOTE: This isn't required to be added immediately in order to proceed with account creation, but is the only way an HCP can view your data within their account.
 - d. Select if you agree to Glooko's Privacy Policy and Terms of Use.
 - e. Click **"Register"**.

- 2 On the next screen, select **"Omnipod DASH Cloud"**, and click **"Next"**.

- 3 Review your sync compatibility, then click **"Next"**.

TIP: You can also create an account from the Glooko Mobile App!

Step 1

***PLEASE NOTE:** if your Healthcare Professional has uploaded your current Omnipod DASH® PDM to Glooko before, or data from another diabetes technology such as CGM or BGM, you may already have a Glooko account. If you can't remember your login credentials, contact Glooko Support to avoid duplicate account creation.


III. Turn On and Connect to Wi-Fi Network

- 1 From the PDM home screen, navigate to the menu icon in the upper left side.
 - a. Go to Settings > PDM Device > Wi-Fi.
 - b. Tap to toggle Wi-Fi on. When the Wi-Fi toggle is set to on, the PDM scans for available Wi-Fi networks.
- 2 To connect to a Wi-Fi network, tap **"Wi-Fi"** to display a list of available Wi-Fi networks. Then tap the name of the desired network and enter the network password.


NOTE: the PDM does not accept networks that require terms and conditions through an Internet browser (*i.e. hotel lobby, Starbucks*). Must be password protected.

 - a. Tap **"Connect"**.
- 3 On the upper right of the screen, tap on the three dots.
 - a. Tap **"Advanced"**.
 - b. Tap **"Keep Wi-Fi on During Sleep"**.
 - c. Tap **"Always"** or **"Only when plugged in"**.
- 4 Check if PDM data is flowing to the Insulet cloud in the **"About"** screen. Check that the **"Last Sync to Cloud Date"** is indeed recent.


If it is NOT, contact Insulet Product Support ASAP.



Step 1



Step 3a



Step 3c

IV. Enable Automatic Upload With Glooko Via Web Browser

- 1 Go to my.glooko.com and sign in to your Glooko account.
- 2 Navigate to the upper right side of the home screen next to your initials, click the down arrow.
 - a. Select **"Settings"**.
 - b. Scroll down to bottom of the page to **"Apps & Devices"**.
 - c. Click **"Connect to DASH Cloud"**.
- 3 A prompt will direct you to Connect to PodderCentral™.
 - a. Click **"Connect"**.
 - b. You will be taken to PodderCentral™ to login.
 - **TIP:** See Section I in this guide for further instructions on how to set up your PodderCentral™ account.
- 4 Once logged into PodderCentral™, you will be prompted to enter the serial number of your new PDM.

TIP: The serial number can be found either on the back of PDM or in the About section of your PDM (Settings --> About).


 - a. Select **"Next"**.
- 5 You will be prompted to review the **"DASH to Glooko Directed Disclosure"**.
 - a. Select **"AGREE"**.
- 6 You will then be returned to Glooko to a screen that states **"Connection Pending"**.
 - a. Click on the checkbox next to **"Yes, share my data with Insulet"** (optional).
 - b. Select **"Done"**.
- 7 Your status will now be updated to **"DASH Cloud: Pending"**.

TIP: To check status in Glooko, return to the **"Settings"** page and scroll down to Apps & Devices to see **"DASH Cloud: Connected"**.


NOTE: If the PDM is new and does not have a complete day of data recorded on it, status will change to **"Connected"** at midnight, as long as it is connected to Wi-Fi. If the PDM has been used with an active Pod for more than a day, then the status will change to **"Connected"** within ~10 minutes, as long as it is connected to Wi-Fi.

Here is a link to a helpful how-to video outlining the steps to automatic upload:


omnipod.com/DASHglooko




Step 2a




Step 2c



Step 3a




Step 3b




Step 4



Step 5




Step 6



Step 7


V. Enable Automatic Upload via Glooko Mobile App

- 1 Open the Glooko mobile app and click **"Sync"** in the top-right.
 - a. Tap **"New Device"**.
 - b. Tap **"Insulin Pumps"**, then tap **"Continue"**.
 - c. Tap **"Omnipod DASH Cloud"**, then tap **"Add"**.
 - d. Tap **"Connect"**.



Step 1

- 2 Follow these on-screen prompts to complete connection:
 - a. Log into PodderCentral account.
 - b. Register PDM.
 - c. Sign the Directed Disclosure Terms and Conditions.
 - d. Navigate back to Glooko.




Step 1a

- 3 The DASH® Cloud account status will show as **"Pending"** while Glooko establishes a connection. Status will change from **"Pending"** to **"Active"** once the connection is established.


NOTE: If the PDM is new and does not have a complete day of data recorded on it, status will change to **"Connected"** at midnight, as long as it is connected to Wi-Fi. If the PDM has been used with an active Pod for more than a day, then the status will change to **"Connected"** within ~10 minutes, as long as it is connected to Wi-Fi.

VI. Connect a ProConnect Code


- 1 Contact your Healthcare Provider to retrieve their unique ProConnect Code.
- 2 Log onto Glooko at my.glooko.com.
- 3 On the upper right side of the home screen next to your initials, select **"Settings"**.
- 4 Scroll down to **"Account"**. Under Primary Provider select **"+ Add New Code"**.
- 5 Select **"Continue"**.
- 6 Enter ProConnect Code.
 - a. Select **"Submit"**.
 - b. Select **"Ok"**.




Step 1




Step 2



Step 3



Step 4



Step 5



For data display, Glooko account creation, Glooko login, or connecting DASH to Glooko status questions, contact Glooko Support at:

1-800-206-6601 • support@glooko.com



For Wi-Fi connection questions, PodderCentral support or PDM/Pod questions, contact Insulet Support at:

1-800-591-3455 • omnipod.com



Insulet Corporation 100 Nagog Park, Acton, MA 01720 • 800-591-3455 • omnipod.com

© 2020 Insulet Corporation. Omnipod, the Omnipod logo, DASH, the DASH logo, Podder and PodderCentral are trademarks or registered trademarks of Insulet Corporation. All rights reserved. Glooko is a trademark of Glooko, Inc. and used with permission. Insulet Corporation does not operate myglooko.com and is not responsible for its conduct or content. Visit glooko.com/privacy/ or call 800.206.6601 for more information. The Bluetooth® word mark and logos are registered trademarks owned by the Bluetooth SIG, Inc. and any use of such marks by Insulet Corporation is under license.
INS-ODS-08-2020-00066 V1.0